



Trauma Catalogue

CANNULATED SCREWS

implants for life™

Table of Contents

3.5 mm Cannulated Screws _____	4
4.0 mm Cannulated Screws _____	5
4.5 mm Cannulated Screws _____	6
Washers _____	7
Cannulated Instruments 3.5-4.0-4.5 _____	8
6.5 mm Cannulated Screws _____	10
7.3 mm Cannulated Screws _____	11
Washers _____	12
Cannulated Instruments 6.5-7.3 _____	13

Screws

3.5mm Cannulated Screws

	Titanium	Stainless Steel	Description
Ref No	41-45035-018	41-25035-018	3.5mm Cannulated Screws, 18mm long, 8mm thread
	41-45035-020	41-25035-020	3.5mm Cannulated Screws, 20mm long, 8mm thread
	41-45035-022	41-25035-022	3.5mm Cannulated Screws, 22mm long, 10mm thread
	41-45035-024	41-25035-024	3.5mm Cannulated Screws, 24mm long, 10mm thread
	41-45035-026	41-25035-026	3.5mm Cannulated Screws, 26mm long, 10mm thread
	41-45035-028	41-25035-028	3.5mm Cannulated Screws, 28mm long, 10mm thread
	41-45035-030	41-25035-030	3.5mm Cannulated Screws, 30mm long, 12mm thread
	41-45035-032	41-25035-032	3.5mm Cannulated Screws, 32mm long, 12mm thread
	41-45035-034	41-25035-034	3.5mm Cannulated Screws, 34mm long, 12mm thread
	41-45035-036	41-25035-036	3.5mm Cannulated Screws, 36mm long, 12mm thread
	41-45035-038	41-25035-038	3.5mm Cannulated Screws, 38mm long, 12mm thread
	41-45035-040	41-25035-040	3.5mm Cannulated Screws, 40mm long, 16mm thread
	41-45035-042	41-25035-042	3.5mm Cannulated Screws, 42mm long, 16mm thread
	41-45035-044	41-25035-044	3.5mm Cannulated Screws, 44mm long, 16mm thread
	41-45035-046	41-25035-046	3.5mm Cannulated Screws, 46mm long, 16mm thread
	41-45035-048	41-25035-048	3.5mm Cannulated Screws, 48mm long, 16mm thread
	41-45035-050	41-25035-050	3.5mm Cannulated Screws, 50mm long, 16mm thread
	41-45035-052	41-25035-052	3.5mm Cannulated Screws, 52mm long, 16mm thread
	41-45035-054	41-25035-054	3.5mm Cannulated Screws, 54mm long, 16mm thread
	41-45035-056	41-25035-056	3.5mm Cannulated Screws, 56mm long, 16mm thread
	41-45035-058	41-25035-058	3.5mm Cannulated Screws, 58mm long, 16mm thread
	41-45035-060	41-25035-060	3.5mm Cannulated Screws, 60mm long, 16mm thread

Description

Cortical cannulated screws partially threaded, indicated for fracture fixation of small bones and small bone fragments, such as metatarsal, arthrodeses of the carpals & phalanges and ligament fixation. The self-drilling and self-tapping design facilitates the ease of insertion. The wide cannulation allows the insertion of thick and stiff guide wires which assist the accurate placement of screw.

Screw Section

*Head diameter 6,0mm
Core diameter 2,5mm
Thread diameter 3,5mm
Hex size 2,5mm
Cannulation 1,2mm*

Material

*Certified Medical Titanium Alloy ISO 5832-3
Certified Medical Stainless Steel ISO 5832-1*



Screws

4.0mm Cannulated Screws

	Titanium	Stainless Steel	Description
Ref No	41-45040-020	41-25040-020	4.0mm Cannulated Screws, 20mm long, 8mm thread
	41-45040-022	41-25040-022	4.0mm Cannulated Screws, 22mm long, 10mm thread
	41-45040-024	41-25040-024	4.0mm Cannulated Screws, 24mm long, 10mm thread
	41-45040-026	41-25040-026	4.0mm Cannulated Screws, 26mm long, 10mm thread
	41-45040-028	41-25040-028	4.0mm Cannulated Screws, 28mm long, 10mm thread
	41-45040-030	41-25040-030	4.0mm Cannulated Screws, 30mm long, 12mm thread
	41-45040-032	41-25040-032	4.0mm Cannulated Screws, 32mm long, 12mm thread
	41-45040-034	41-25040-034	4.0mm Cannulated Screws, 34mm long, 12mm thread
	41-45040-036	41-25040-036	4.0mm Cannulated Screws, 36mm long, 12mm thread
	41-45040-038	41-25040-038	4.0mm Cannulated Screws, 38mm long, 12mm thread
	41-45040-040	41-25040-040	4.0mm Cannulated Screws, 40mm long, 16mm thread
	41-45040-042	41-25040-042	4.0mm Cannulated Screws, 42mm long, 16mm thread
	41-45040-044	41-25040-044	4.0mm Cannulated Screws, 44mm long, 16mm thread
	41-45040-046	41-25040-046	4.0mm Cannulated Screws, 46mm long, 16mm thread
	41-45040-048	41-25040-048	4.0mm Cannulated Screws, 48mm long, 16mm thread
	41-45040-050	41-25040-050	4.0mm Cannulated Screws, 50mm long, 16mm thread
	41-45040-052	41-25040-052	4.0mm Cannulated Screws, 52mm long, 16mm thread
	41-45040-054	41-25040-054	4.0mm Cannulated Screws, 54mm long, 16mm thread
	41-45040-056	41-25040-056	4.0mm Cannulated Screws, 56mm long, 16mm thread
	41-45040-058	41-25040-058	4.0mm Cannulated Screws, 58mm long, 16mm thread
	41-45040-060	41-25040-060	4.0mm Cannulated Screws, 60mm long, 16mm thread

Description

Cancellous cannulated screws partially threaded, indicated for fracture fixation of small bones and small bone fragments. The self-drilling and self-tapping design facilitates the ease of insertion. The wide cannulation allows the insertion of thick and stiff guide wires which assist the accurate placement of screw.

Screw Section

*Head diameter 6,0mm
Core diameter 2,5mm
Thread diameter 4,0mm
Hex size 2,5mm
Cannulation 1,2mm*

Material

*Certified Medical Titanium Alloy ISO 5832-3
Certified Medical Stainless Steel ISO 5832-1*



Screws

4.5mm Cannulated Screws

	Titanium	Stainless Steel	Description
Ref No	41-45045-030	41-25045-030	4.5mm Cannulated Screws, 30mm long, 16mm thread
	41-45045-032	41-25045-032	4.5mm Cannulated Screws, 32mm long, 16mm thread
	41-45045-034	41-25045-034	4.5mm Cannulated Screws, 34mm long, 16mm thread
	41-45045-036	41-25045-036	4.5mm Cannulated Screws, 36mm long, 16mm thread
	41-45045-038	41-25045-038	4.5mm Cannulated Screws, 38mm long, 16mm thread
	41-45045-040	41-25045-040	4.5mm Cannulated Screws, 40mm long, 16mm thread
	41-45045-042	41-25045-042	4.5mm Cannulated Screws, 42mm long, 16mm thread
	41-45045-044	41-25045-044	4.5mm Cannulated Screws, 44mm long, 16mm thread
	41-45045-046	41-25045-046	4.5mm Cannulated Screws, 46mm long, 16mm thread
	41-45045-048	41-25045-048	4.5mm Cannulated Screws, 48mm long, 16mm thread
	41-45045-050	41-25045-050	4.5mm Cannulated Screws, 50mm long, 16mm thread
	41-45045-052	41-25045-052	4.5mm Cannulated Screws, 52mm long, 16mm thread
	41-45045-054	41-25045-054	4.5mm Cannulated Screws, 54mm long, 16mm thread
	41-45045-056	41-25045-056	4.5mm Cannulated Screws, 56mm long, 16mm thread
	41-45045-058	41-25045-058	4.5mm Cannulated Screws, 58mm long, 16mm thread
	41-45045-060	41-25045-060	4.5mm Cannulated Screws, 60mm long, 16mm thread
	41-45045-062	41-25045-062	4.5mm Cannulated Screws, 62mm long, 16mm thread
	41-45045-064	41-25045-064	4.5mm Cannulated Screws, 64mm long, 16mm thread

Description

Cortical cannulated screws partially threaded, indicated for fracture fixation of large bones and large bone fragments. The self-drilling and self-tapping design facilitates the ease of insertion. The wide cannulation allows the insertion of thick and stiff guide wires which assist the accurate placement of screw.

Screw Section

*Head diameter 6,0mm
Core diameter 3,0mm
Thread diameter 4,5mm
Hex size 2,5mm
Cannulation 1,6mm*

Material

*Certified Medical Titanium Alloy ISO 5832-3
Certified Medical Stainless Steel ISO 5832-1*



Washers

Washer

	<i>Titanium</i>	<i>Stainless Steel</i>	<i>Description</i>
<i>Ref No</i>	60-45201-010	60-25201-010	Washer Ø3.5mm
	60-45201-020	60-25201-020	Washer Ø4.0mm
	60-45201-030	60-25201-030	Washer Ø4.5mm



Instruments



VA-58601605
Cannulated Drill Bit ø 3.0

VA-58501630
Drill Bit Stop



VA-58601628
Screw Extractor



VA-58601614
Cannulated Tap



VA-58601601
Guide Pin ø 1.2



VA-58601624
Cleaning Guide Pin



VA-58601609
Cannulated Countersink Drill



VA-58601618
Cannulated Tap ø 4.0



VA-58601607
Cannulated Drill Bit ø 2.5



VA-58601616
Cannulated Tap ø 4.5



VA-58691612
Drill Sleeve ø 2.5



VA-58601613
Guide Pin Sleeve ø 1.2



VA-58601637
Stop Wrench



VA-58601635
Guide Pin Depth Gauge



VA-58601611
Drill Sleeve ø 3.0



VA-58601620
Guide Pin Sleeve

VA-58601012
Screw Holding Forceps




VA-58601632
Cannulated Screwdriver SW2.5



VA-58601027
Hex Screwdriver SW2.5

VA-58601622
Screw Holding Sleeve



Screws

6.5mm Cannulated Screws

	Titanium	Stainless Steel	Description
Ref No	41-45065-050	41-25065-050	6.5mm Cannulated Screws, 50mm long, 18mm thread
	41-45065-055	41-25065-055	6.5mm Cannulated Screws, 55mm long, 18mm thread
	41-45065-060	41-25065-060	6.5mm Cannulated Screws, 60mm long, 18mm thread
	41-45065-065	41-25065-065	6.5mm Cannulated Screws, 65mm long, 18mm thread
	41-45065-070	41-25065-070	6.5mm Cannulated Screws, 70mm long, 18mm thread
	41-45065-075	41-25065-075	6.5mm Cannulated Screws, 75mm long, 18mm thread
	41-45065-080	41-25065-080	6.5mm Cannulated Screws, 80mm long, 18mm thread
	41-45065-085	41-25065-085	6.5mm Cannulated Screws, 85mm long, 18mm thread
	41-45065-090	41-25065-090	6.5mm Cannulated Screws, 90mm long, 18mm thread
	41-45065-095	41-25065-095	6.5mm Cannulated Screws, 95mm long, 18mm thread
	41-45065-100	41-25065-100	6.5mm Cannulated Screws, 100mm long, 18mm thread
	41-45065-105	41-25065-105	6.5mm Cannulated Screws, 105mm long, 18mm thread
	41-45065-110	41-25065-110	6.5mm Cannulated Screws, 110mm long, 18mm thread

Description

Cancellous cannulated screws partially threaded, indicated for fracture fixation of large bone and large bone fragments, it is particularly intended for femoral neck fracture, slipped capital femoral epiphysis, tibial plateau fracture, ankle arthodeses, intercondylar femur fracture and sacroiliac joint disruption. The self-drilling and self-tapping design facilitates the ease of insertion. The wide cannulation allows the insertion of thick and stiff guide wires which assist the accurate placement of screw.

Screw Section

*Thread Diameter 6.5
Cannulation 2.6
Hex Size 4.0
Head Diameter 8.0*

Material

*Certified Medical Titanium Alloy ISO 5832-3
Certified Medical Stainless Steel ISO 5832-1*



Screws

7.3mm Cannulated Screws

	Titanium	Stainless Steel	Description
Ref No	42-45073-050	42-25073-050	7.3mm Cannulated Screws, 50mm long, 18mm thread
	42-45073-055	42-25073-055	7.3mm Cannulated Screws, 55mm long, 18mm thread
	42-45073-060	42-25073-060	7.3mm Cannulated Screws, 60mm long, 18mm thread
	42-45073-065	42-25073-065	7.3mm Cannulated Screws, 65mm long, 18mm thread
	42-45073-070	42-25073-070	7.3mm Cannulated Screws, 70mm long, 18mm thread
	42-45073-075	42-25073-075	7.3mm Cannulated Screws, 75mm long, 18mm thread
	42-45073-080	42-25073-080	7.3mm Cannulated Screws, 80mm long, 18mm thread
	42-45073-085	42-25073-085	7.3mm Cannulated Screws, 85mm long, 18mm thread
	42-45073-090	42-25073-090	7.3mm Cannulated Screws, 90mm long, 18mm thread
	42-45073-095	42-25073-095	7.3mm Cannulated Screws, 95mm long, 18mm thread
	42-45073-100	42-25073-100	7.3mm Cannulated Screws, 100mm long, 18mm thread
	42-45073-105	42-25073-105	7.3mm Cannulated Screws, 105mm long, 18mm thread
	42-45073-110	42-25073-110	7.3mm Cannulated Screws, 110mm long, 18mm thread
	42-45073-115	42-25073-115	7.3mm Cannulated Screws, 115mm long, 18mm thread

Description

Cancellous cannulated screws partially threaded, indicated for fracture fixation of large bone and large bone fragments, it is particularly intended for femoral neck fracture, slipped capital femoral epiphysis, tibial plateau fracture, ankle arthodeses, intercondylar femur fracture and sacroiliac joint disruption. The self-drilling and self-tapping design facilitates the ease of insertion. The wide cannulation allows the insertion of thick and stiff guide wires which assist the accurate placement of screw.

Screw Section

Thread Diameter 7.3

Cannulation 2.6

Hex Size 4.0

Head Diameter 8.0

Material

Certified Medical Titanium Alloy ISO 5832-3

Certified Medical Stainless Steel ISO 5832-1



WASHERS

Washer

	<i>Titanium</i>	<i>Stainless Steel</i>	<i>Description</i>
<i>Ref No</i>	60-45201-040	60-25201-040	Washer Ø6.5mm
	60-45201-050	60-25201-050	Washer Ø7.3mm



Instruments



VA-58601718
Cleaning Guide Pin



VA-58601734
Guide Pin Depth Gauge



VA-58601703
Cannulated Drill Bit



VA-58601718
Cannulated Tap ø 7.3

VA-58601725
Drill Bit Stop



VA-58601719
Cannulated Tap ø 6.5



VA-58601736
Guide wrench



VA-58601701
Guide Pin



VA-58601732
Hex Screwdriver SW4.0



VA-58601731
Cannulated Hex Screwdriver



VA-58601729
Multi-Holes Guide



VA-58601714
Two Hole Guide



VA-58601727
Screw Extractor



VA-58601705
Cannulated Countersink Drill




VA-58601707
Drill Sleeve


VA-58601716
Guide Pin Sleeve



VA-58601711
Support Device

58601709
Dilator





1 Kings Avenue, London N21 3 NA, UK
t: +44 (0) 208 3707764 f: +44 (0) 208 3707292
website: www.vadinimplants.com
e-mail: info@vadinimplants.com